

Span F500S

Powered Sit-to-Stand Patient Lift



For users
Requiring
**MODERATE
ASSISTANCE**

- **High Capacity:** 500 lb (228 kg) safe working load.
- **Care Flexibility:** Underbed clearance allows use with most low beds.
- **Safe and Reliable:** Redundant controls; overload protection; emergency lowering; 3-year frame warranty.
- **Easy operation:** Simple charging; audible/visual indicators; locking rear casters.

SPAN
better science. better care.



LIFTS



SURFACES



BEDS



SEATING



FALL PROTECTION



POSITIONERS



SKIN CARE



F500S Sit-To-Stand Patient Lift

An essential tool for safe, efficient patient transfers between bed, chair and commode. Compact and easy to maneuver in a variety of settings, the F500S delivers a winning combination of safety, simplicity, and value.



LCD battery level display; audible low-battery indicator



2-button pendant with LED indicator



Underbed clearance, allows use with most low beds



Contoured knee pad with 3 height settings

| SPECIFICATIONS | |
|---------------------------------------|------------------|
| Weight Capacity | 500 lb / 228 kg |
| External base width (closed) | 24.9" |
| Internal leg width (opened) | 33.9" |
| Upright mast height | 40" |
| Overall base length | 37.5" |
| Base Height | 4.9" |
| Min. lifting height with spreader bar | 38" |
| Max. lifting height with spreader bar | 60.5" |
| Knee pad length extension range | 4.5" |
| Max. boom width | 17.7" |
| Net Weight | 111 lb / 50.3 kg |
| Shipping Weight | 119 lb / 53.9 kg |

- **Easy base width adjustment** with locked position; accommodates chairs and commodes.
- **Wide non-slip footplate** provides stable base for transfer.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Backlit LCD display** during operation for ease of viewing.
- **Simple, fast assembly** using tools provided.
- **Certified to:** ISO 10535; IEC/EN/ANSI/AAMI/CAN/CSA-C22.2 60601-1, 60101-1-2, and 60601-1-11; SGS Mark for Safety.

| Item # | Description |
|------------------|--|
| F500S | Sit-To-Stand Patient Lift |
| Warranty: | 3 years (frame and spreader bar); 2 years (electronics); 1 year (battery) |

Part of Span's line of Safe Patient Handling and Mobility solutions, including:

F500P Full Body Patient Lift



Floor Lifts

F450T Standing Transfer Aid



Fixed Lithium-ion



Overhead Lifts

Portable Lithium-ion

